



```
% design an appropriate filter for ECG signal  
% detection
```

```
% step 1: create noise ECG signal  
load noisyecg.mat  
fs = 360;
```

```
% step 2: plot ECG signal  
plot (noisyECG_withTrend);  
title("NOISY ECG SIGNAL");  
xlabel("Time");  
ylabel("Amplitude");
```



```
% step 3: design a 60hz eliminator (low pass filter) using the filter  
% taskblock from TASK option  
% Design a digital filter  
lowpass1 = designfilt('lowpassfir', ...  
    'FilterOrder',20,'CutoffFrequency',70, ...  
    'PassbandRipple',1,'StopbandAttenuation',60, ...  
    'SampleRate',fs);  
  
% Visualize magnitude and phase responses  
freqz(lowpass1.Coefficients,1,[],fs)
```




% step 4: apply the noisy ECG singal to design singal to the design filter and the save enhanced singnal

```
Enhanced_singal = filter(lowpass1,noisyECG_withTrend);
```

% step 5: plot the enhanced singnal & observe the singal give it a title

```
plot(Enhanced_singal);
title("Enhanced ECG singal");
xlabel("Time");
ylabel("Amplitude")
```



% step 6: now the bandpass filter for the range of 0.25 to 40hz

% Design a digital filter

```
band2 = designfilt('bandpassfir', ...
```

```

'FilterOrder',30,'CutoffFrequency1',0.25, ...
'CutoffFrequency2',40,'StopbandAttenuation1',70, ...
'PassbandRipple',1,'StopbandAttenuation2',70, ...
'SampleRate',fs);

```

```


% Visualize magnitude and phase responses

```

```

freqz(band2.Coefficients,1,[],fs)

```



```

% step 7: apply the designed filter to enhanced singal
filtered_signal = filter(band2, Enhanced_singal);

```

```

% step 8: plot the filter singal and observe the result
plot(filtered_signal);
title("Filtered ECG singal");
xlabel("Time");
ylabel("Amplitude")

```

